



UVBIOTEK™ MOBILE-LITE™
6 NARROWBAND UVB LAMPS

PATENT PENDING

The NEW Mobile-Lite™: Phototherapy right where you need it, when you need it

Right for your lifestyle

The Mobile-lite™ is the first completely portable phototherapy system designed specifically to treat hands and/or feet. It collapses down to the size of a large briefcase, has a built-in handle, and only weighs about 40 lbs, so it can be carried wherever you go.



Right for your hands

There is no safer or more effective way to treat psoriasis build-up on the palms of the hands. The Mobile-Lite™ can quickly be set up on any flat surface to provide the necessary treatment for the largest number of sufferers.



Right for your feet

If you have psoriasis flaring on the soles of your feet, you know how painful standing or walking can be. The more you can stay off your feet, the sooner you will be able to get clear again. The Mobile-Lite™ sets up just as easily on the floor to treat the feet while you keep your seat.



At home, at work, on vacation—UVB therapy that's there for you.

For nearly 100 years, UVB phototherapy has been in medical use as one of the safest and most effective ways to treat psoriasis. Unlike medications which treat the condition from the inside out, phototherapy treats the skin directly from the outside, and with virtually no side effects.

Psoriasis affecting the palms and the soles of the feet can be very difficult to treat using conventional phototherapy systems. This is because very little UVB can penetrate through the thick plaques that build-up on the hands, and the bottoms of the feet get no UVB at all when you are standing up.

For this reason, UVBioTek™ developed the Mobile-Lite™ system to focus the therapeutic light directly where it is needed most. And, because it can be used just about anywhere and at anytime, it is easy to follow the recommended treatment schedule of every other day.

Full-body phototherapy can take from 6 months to a year to clear up hand or foot plaques. Our compact, portable Mobile-Lite™ can usually get you clear again in a fraction of that time. The Mobile-Lite™ does your hands and feet right.

We offer a variety of affordable financing plans to fit your budget, and can help you with getting insurance approval. Call 800.UVBioTek (800.882.4683).

 **UVBIOTEK™**
PHOTOTHERAPY DESIGNED TO FIT™

The NEW Mobile-Lite™

Flip-top lid closes to help protect the lamps while transporting, and helps contain the UVB light while treating.



Remote control unit

6 narrowband lamps with acrylic shield.

Power cord storage area.

System shown with the vinyl UV shield removed.
Patent Pending

Wherever you go, your treatment can come with you.

You just place the Mobile-Lite™ on a desk, a table, or the floor, open the case, and pull the vinyl UV shield forward and secure it into position.

Our unique shield minimizes the amount of UV light that escapes into your room.

With a timer that can be used outside the UV shield and a 6-foot power cord built right in — you're ready to go.

The UVBioTek™ Mobile-Lite™ can also be turned on its side and used, with or without the UV shield, to treat almost any small area of the body. So you can get focused, effective phototherapy wherever and whenever you need it.



The right fit for them

“Results are spectacular. Completely clear... Everyone with psoriasis should own one.” —M.A., Massachusetts

“My treatment is progressing well. I’ve gone from full blown flare up to nearly complete remission. Awesome results... I could not be more pleased... cannot wait for my dermatologist to see the results...” —S.M., Colorado

Home system specifications

| | |
|--|---|
| UVB Lamps | 6 Narrowband UVB lamps |
| Dimensions (closed) Dimensions (open) | 5" x 34 1/2" x 17 1/2" (HWD) 18" x 34 1/2" x 17 1/2" (HWD) |
| Power Needs | 120VAC @ 60Hz (std US household power); minimum 15amp breaker |
| Power Cord | 6 foot coiled power cord / standard three-prong plug |
| Shipping Weight | 50 lbs. |

Get a Mobile-Lite™ system and start getting phototherapy treatments where and when you need them.

To see if UVBioTek™ is the right fit for you, call **800.UVBioTek** (800.882.4683) or visit www.uvbiotek.com today!

CAUTION: Federal law restricts this device for sale by or on the order of a physician.

