


UVBIOTEK™  
 MODEL 100B  
 10 NARROWBAND  
 UVB LAMPS  
 PATENT PENDING

# The NEW Model 100B: Phototherapy that offers you more for less.

## More space and style



Unit closes to conceal the lamps when not in use. And door panels come in your choice of attractive wood grains and solid colors that will blend in perfectly with the style of any room.

## More therapeutic benefits



With the 100B light panel, you treat both sides of the body separately—first the front, then the back. The door panels swing out to reveal custom reflector panels that redirect more UVB light to the sides of the body. This can help even out the exposure.

## More for your money



While treating one side of the body at a time means longer sessions, the 100B produces quality therapeutic results in a compact, affordably priced system.

## Now you can enjoy all the benefits of home phototherapy, without sacrificing your décor—or your wallet.

The medicinal history of UVB phototherapy dates back almost 100 years. In that time, it has been found generally safe and effective in treating a wide range of conditions like psoriasis, eczema, and pruritis. It is effective for the majority of patients, and doesn't carry the long list of side effects usually encountered when taking medications. More recently, narrowband UVB lamps have made phototherapy even more effective for a larger number of people, and for a larger number of conditions. For vitiligo,


in particular, it is fast becoming the only truly effective treatment available.

Now, with the Model 100B single panel system, phototherapy treatments in the convenience of your own home are easy and more affordable than ever. Home phototherapy gives you a flexible schedule which allows you to get your treatment every other day, as recommended. And, since home systems are covered by most insurance, it may be even easier than you think to get one today.

We offer a variety of affordable financing plans to fit your budget, and can help you with getting insurance approval. Call 800.UVBioTek (800.882.4683).




# The NEW Model 100B




10 narrowband lamps with acrylic shield.

Handsome doors conceal the system from view; you choose the look you prefer.



Mounts to virtually any wall; plugs into any wall outlet.

Door panels have a reflective interior to redirect more UVB light to the body sides.



Patent Pending

## The right fit for them

"I am EXTREMELY impressed with your product and its results so far... The narrowband bulbs have made an incredible difference... I am truly blown away with how quickly I have gotten results and am kicking myself for not purchasing this unit years ago." — R.F., California

"...everything is great. The system works beautifully. There is a big improvement in my condition. I think I should be clear in another few weeks... Your company did a great job designing this cabinet for looks and home use." — L.C., New Jersey

## Home system specifications

UVB Lamps	10 Narrowband lamps
Dimensions	86" x 22" x 5" (HWD)
Power Needs	120VAC @ 60Hz (std US household power); minimum 15amp breaker
Power Cord	10 foot cord with standard three-prong plug
Shipping Weight	180 lbs.

With a design that will fit in your home and a price you can afford, the UVBioTek™ Model 100B is a bright idea for you.

To see if UVBioTek™ is the right fit for you, call **800.UVBioTek** (800.882.4683) or visit [www.uvbiotek.com](http://www.uvbiotek.com) today!

**CAUTION:** Federal law restricts this device for sale by or on the order of a physician.

